




# Rocket® Fetal Blood Sampling Kits

Rocket® FBS Kits offer simple, safe and cost effective solutions for routine fetal blood sampling.

- **RELIABLE SAMPLING:** In conjunction with the world's leading manufacturers of blood gas analysers, Rocket Medical have specifically developed kits compatible with the following systems:
  - ✓ Siemens® (Bayer) RapidLab™ Series
  - ✓ Radiometer® ABL™ 800, 700, 77, and 5/500 Series
  - ✓ Instrumentation Laboratory® GEM™
  - ✓ ROCHE micro volume systems.
- **COMPREHENSIVE KITS:** each designed to complete up to 5 individual FBS samplings over the delivery period. Each kit is ready for immediate use to minimise delay when time is precious.
- **AMNILUME™ COMPATIBLE:** For maximum illumination with minimum of fuss, all Rocket® FBS Kits are designed for use with the revolutionary **Rocket® Amnilume™ Single Use Light Source**
- **REDUCE RISK AND STRESS** caused by incompatible sampling tubes and missing or worn instruments frequently found in reusable trays.
- **EFFICIENT AND COST EFFECTIVE** when compared to reprocessing and replacing worn or lost reusable instruments.



- ✓ Endoscope & Obturator
  - ✓ Amnilume™ Single Use Light Source
  - ✓ 6 x 12" cotton swab sticks
  - ✓ 5 x gauze swabs
  - ✓ 1 x disposable sponge holder
  - ✓ 2 x Rocket® FBS Blades
  - ✓ 1 x Rocket® FBS Blade Holder
  - ✓ 1 x Yellow Soft Paraffin BP
  - ✓ 1 x capillary tube holder
  - ✓ 10 x capillary tube stoppers
  - ✓ Tray wrapped in waterproof drape
- Packed individually in cartons of 10 units

Analyser	Recommended volume	Rocket® FBS Kit	Rocket® FBS Kit with Amnilume™
SIEMENS® (Bayer) Rapidlab™ Series	100µl	-	R57026
ABL™ 800 FLEX with FLEXMODE selected	>35µl	Any	R57012-LUME R57018-LUME R57024-LUME
ABL™ 700 Series with tube volume selected	>35µl	R57018	R57018-LUME
ABL™ 77 Series	85µl	R57024	R57024-LUME
ABL™ 5/500 Series	85µl	R57024	R57024-LUME
Instrumentation Laboratory (IL) GEM™ Series	85µl	R57024	R57024-LUME
ROCHE® Systems	55µl	R57018	R57018-LUME

<b>In the UK - call Customer Services: 0191 419 6988</b>	<b>Outside of UK - call Customer Services: +44 192 365 1404</b>
--	---